



Product List 2017



Luminex System

Luminex

<i>Luminex</i>	LX200 IS	Luminex 200 IS - System	Analyzer, XY-Platform, Sheath Delivery System, PC, xPonent 3.1 Software, Barcodescanner. Detect up to 100 individual analytes. Compatible with all listed assays
----------------	----------	--------------------------------	--

Optiplex Assays

Borrelia

Luminex

<i>Diamex</i>	IN0503 201004	Borrelia IgM	Fluorescent Bead Assay (FBA) for the detection of IgM against a Borrelia lysate and 8 recombinant Borrelia afzelii/garinii antigens (OspA, OspC (2x), p100, p18, p39, p41-I, VlsE-C6)	96 tests
<i>Diamex</i>	IN0504 201003	Borrelia IgG	Fluorescent Bead Assay (FBA) for the detection of IgG against a Borrelia lysate and 8 recombinant Borrelia afzelii/garinii antigens (OspA, OspC, p100, p18, p39, p41-I, p58, VlsE-C6)	96 tests
<i>Diamex</i>	IN0501	Borrelia IgM Screen	Fluorescent Bead Assay (FBA) for the detection of IgM against a lysate of Borrelia afzelii and recombinant Osp C of Borrelia afzelii and garinii	96 tests
<i>Diamex</i>	IN0502	Borrelia IgG Screen	Fluorescent Bead Assay (FBA) for the detection of IgG against 2 lysate antigens (Borrelia afzelii and garinii)	96 tests

HPV

Luminex

<i>Diamex</i>	IN0601	HPV Genotyping	Fluorescent Bead Assay (FBA) for the identification of 24 high and low-risk genital HPV genotypes in PCR amplified samples	96 tests
---------------	--------	-----------------------	--	----------

FIDIS Assays

Autoimmune Diagnostics				
Luminex				
Theradiag	MX 002 411002	Thyro (TPO, TG)	Simultaneous detection of autoantibodies against Thyroperoxidase (TPO) and Thyroglobulin (TG)	96 tests
Theradiag	MX 005 411005	dsDNA	Quantitative detection of autoantibodies directed against double stranded DNA (dsDNA)	96 tests
Theradiag	MX 006 411006	Connective 10	Simultaneous detection of autoantibodies against dsDNA, SS-A (60 kDa), TRIM21 (previously named SS-A 52 kDa), SS-B, Sm, Sm/RNP, Scl-70, Jo-1, Ribosomal Proteins, Centromere (CENP-B)	96 tests
Theradiag	MX 117 411117	Connective Profile	Simultaneous detection of autoantibodies against dsDNA, SS-A (60 kDa), TRIM21 (previously named SS-A 52 kDa), SS-B, Sm, Sm/RNP, Scl-70, Jo-1, Ribosomal Proteins, Centromere (CENP-B), PmSCI, PCNA, Histone, U1RNP	96 tests
Theradiag	MX 024	Hep2 Screen	Qualitative detection of anti-nuclear antibodies (IgG) against components of Hep2 cells, ENA's and other cytoplasmic antigens	96 tests
Theradiag	MX 007 411007	Vasculitis (MPO, PR3, GBM)	Simultaneous detection of autoantibodies against Myeloperoxidase (MPO), Proteinase 3 (PR3) and Glomerular Basement Membrane (GBM)	96 tests
Theradiag	MX 012 411012	APS IgG	Simultaneous detection of IgG autoantibodies against Cardiolipin, beta2-Glycoprotein I (β 2GPI), Prothrombin also called Factor II (FII)	96 tests
Theradiag	MX 013 411013	APS IgM	Simultaneous detection of IgM autoantibodies against Cardiolipin, beta2-Glycoprotein I (β 2GPI), Prothrombin also called Factor II (FII)	96 tests
Theradiag	MX 030-48 411030	Cyto Panel	Simultaneous detection of autoantibodies against ATPase H ⁺ /K ⁺ of the Gastric Parietal Cells (PCA), Intrinsic Factor (IF), Ribosomal Proteins, Mitochondrial M2 Proteins, LKM1 and LC1	48 tests
Theradiag	MX 103 411103	Rheuma (RF human and animal)	Simultaneous detection of IgM class Rheumatoid Factor directed against Fc determinants from human or animal IgG Immunoglobulins.	96 tests
Theradiag	MX 104 411104	Anti-CCP	Simultaneous detection of autoantibodies directed against 4 Synthetic Citrullinated Peptides.	96 tests
Theradiag	MX 109 411109	Celiac DPG IgG	Simultaneous detection of IgG autoantibodies against Gliadin, Deamidated Peptides of Gliadin (DPG) and Tissue Transglutaminase (tTG)	96 tests
Theradiag	MX 110 411110	Celiac DPG IgA	Simultaneous detection of IgA autoantibodies against Gliadin, Deamidated Peptides of Gliadin (DPG) and Tissue Transglutaminase (tTG)	96 tests

AtheNA Multi-Lyte Assays

EBV				
Luminex				
<i>Zeus Scientific</i>	A92101G 201032	EBV IgG Plus	Multiplex test for the detection of IgG against VCA, EBNA and EA	96 tests
<i>Zeus Scientific</i>	A92201M 201031	EBV IgM Plus	Multiplex test for the detection of IgM against VCA	96 tests
HSV				
Luminex				
<i>Zeus Scientific</i>	A97101G	HSV 1&2 IgG Plus	Multiplex test for the detection of IgG against HSV-1 and HSV-2 (recombinant gG-1 and gG-2)	96 tests
Measles, Mumps, Rubella, Varizella				
Luminex				
<i>Zeus Scientific</i>	A93111G	MMV IgG Plus	Multiplex test for the detection of IgG against Measles virus, Mumps virus and VZV	96 tests
<i>Zeus Scientific</i>	A93101G	MMRV IgG Plus	Multiplex test for the detection of IgG against Measles, Mumps, Rubella virus and VZV	96 tests
Toxoplasma, Rubella, CMV, HSV				
Luminex				
<i>Zeus Scientific</i>	A81101G 201018	ToRCH IgG Plus	Multiplex test for the detection of IgG against Toxoplasma, Rubella, CMV, HSV-1 and HSV-2	96 tests

Quality Control Sera

Luminex

<i>Theradiag</i>	HM060-2X05 412060	Connective 1	dsDNA, SS-A 60kDa, TRIM 21, SS-B, Sm, Sm/RNP & U1-RNP	2x 0.5 ml
<i>Theradiag</i>	HM061-2X05 412061	Connective 2	dsDNA, Scl-70, Centromere protein B, Jo-1 & Ribosomal proteins	2x 0.5 ml
<i>Theradiag</i>	HM062-2X05 412062	APS IgG	Cardiolipin IgG, β 2 GPI IgG, Prothrombin IgG & Phospholipids IgG	2x 0.5 ml
<i>Theradiag</i>	HM063-2X05 412063	APS IgM	Cardiolipin IgM, β 2 GPI IgM, Prothrombin IgM & Phospholipids IgM	2x 0.5 ml
<i>Theradiag</i>	HM064-2X05 412064	Rheuma	Animal and Human RF IgM & CCP IgG	2x 0.5 ml
<i>Theradiag</i>	HM065-2X05 412065	Thyro	TPO, TG & R-TSH	2x 0.5 ml
<i>Theradiag</i>	HM066-2X05 412066	Cyto	Ribosomal proteins, Mitochondria, LKM-1 & LC-1	2x 0.5 ml
<i>Theradiag</i>	HM067-2X05 412067	Vasculitis	dsDNA, MPO, PR3 & GBM	2x 0.5 ml
<i>Theradiag</i>	HM068-2X05 412068	Gastritis	Gastric Parietal Cells (PCA) & Intrinsic Factor (IF) antibodies	2x 0.5 ml
<i>Theradiag</i>	HM069-2X05 412069	Celiac G	Transglutaminase IgG, Gliadin IgG & Deamidated Gliadin IgG	2x 0.5 ml
<i>Theradiag</i>	HM070-2X05 412070	Celiac A	Transglutaminase IgA, Gliadin IgA & Deamidated Gliadin IgA	2x 0.5 ml
<i>Theradiag</i>	HM071-02X05 412071	ASCA G + A	Anti-Saccharomyces Cerevisiae Antibodies IgA & IgG	2x 0.5 ml
<i>Theradiag</i>		Info:	all ImmunoTrol quality controls are also available in quantities of 100x 0.5 ml. ImmunoTrol controls are all CE marked, ready to use with an excellent stability and suitable for DOT, ELISA, FIDIS (Luminex) and other techniques.	

Accessories for Luminex System

Luminex Accessories				
Luminex				
Luminex	MM1101 201081	Sheath Fluid 20L		20 L
Luminex	MM1112 201080	Sheath Fluid concentrate (20x)		1 L
Luminex	MM1102 201082	Calibrator Cal1	Luminex 100/200 Classification Calibration Microspheres	5 ml
Luminex	MM1103 201083	Calibrator Cal2	Luminex 100/200 Reporter Calibration Microspheres	5 ml
Luminex	MM1104 201084	Control Con1	Luminex 100/200 Classification Control Microspheres	5 ml
Luminex	MM1105 201085	Control Con2	Luminex 100/200 Reporter Control Microspheres	5 ml
Luminex	MM1120 201078	Lx100/200 Calibration Kit		25 uses
Luminex	MM1121 201079	Lx100/200 Performance Verification Kit		25 uses
Zeus Scientific	101D 368001	AtheNA Multi-Lyte 96 Well Filterplate		1 pce
Zeus Scientific	AX9001 201094	AtheNA Multi-Lyte Daily Maint. Reagents	provides all of the materials necessary to perform 24 daily startup and 24 daily shutdown procedures on the Luminex system	
Zeus Scientific	005M 201098	Sample Diluent	for EBV samples with INV:NSC	30 ml
Diamex	MK0101 201090	Optiplex LX Cleaning Kit	All materials and solutions needed for 24 daily startup and shutdown procedures of Luminex LX100/LX200 Analyzers and 2x sample probe needle cleansing procedures	1 pce
Diamex	MM1111 201091	Microtiterplate C8 Polysorb		1 pce
Theradiag	MXP051 411190	Microplate with filtering membrane 0.22µm		1 pce
Theradiag	MXP051 411191	Microplate with filtering membrane 0.22µm		5 pce

Asuragen PCR Tests

C9orf72 Gene

Molecular Biology

Asuragen	49581 380140	AmplideX PCR/CE C9orf72	Reliable and reproducible method for high-resolution genotyping of GGGGCC hexanucleotide repeat expansions in the C9orf72 gene, which is of increasing interest in frontotemporal dementia (FTD) and amyotrophic lateral sclerosis (ALS)	50 rxn
----------	-----------------	--------------------------------	--	--------

Fragile X

Molecular Biology

Asuragen	76008 380103	AmplideX PCR/CE FMR1	Reliably analyze CGG trinucleotide repeats and detect interrupting AGG sequences in the FMR1 gene	100 rxn
Asuragen	49442 380110	AmplideX FMR1 mPCR	Innovative PCR-only approach for the detection of allele-specific methylation in the FMR1 gene and elimination of the need for Southern Blot.	24 rxn
Asuragen	49513 380120	AmplideX FMR1 PCR Process Control	Mixture of genomic DNA containing multiple alleles. Can be used as process or batch control covering multiple size ranges in the FMR1 gene. Well characterized and ready to use (24 µl).	12 rxn
Asuragen	49514 380121	AmplideX FMR1 Methylation & Sensitivity Control	Well characterized, ready to use (24 µl) genomic DNA control. Covers multiple size ranges and different methylation states in the FMR1 gene.	3 rxn

Leukemia

Molecular Biology

Asuragen	86003	Quantidex qPCR BCR-ABL IS	Quantitative real time PCR based assay for monitoring the BCR-ABL fusion gene levels of patients with Chronic Myelogenous Leukemia (CML). All the reagents needed for RT and qPCR step are included. The Kit is a sensitive and multiplexed assay that provides a simpler workflow, direct IS alignment and robust performance especially at low % IS ratios.	60 rxn
----------	-------	----------------------------------	---	--------

NGS (Next Generation Sequencing)

Molecular Biology

Asuragen	49539	Quantidex DNA Assay	Multiplexed quantitative PCR assay that measures the absolute copy number of PCR-amplifiable DNA in a sample and reports PCR inhibition. The Quantidex Assay determines the functional quality of sample DNA using the QFI™ Score, which is the fraction of total genomic DNA copies that can be amplified.	100 rxn
Asuragen	49560	Quantidex NGS Pan Cancer Kit	Multiplexed profiling of the most clinically relevant mutations in a wide range of human cancers with high performance in low quality/quantity samples.	48 rxn

Molecular Biology

Meridian Bioscience		illumigene	<i>illumigene is a simplified new technology platform, based upon Loop Mediated Isothermal Amplification (LAMP) that makes molecular biology testing possible for any laboratory.</i> <i>The very compact illumipro-10 incubator/reader can run up to 10 specimens and results are available in less than</i>	
Meridian Bioscience	280050 308520	illumigene C. difficile	DNA amplification assay for the detection of toxigenic Clostridium difficile	50 tests
Meridian Bioscience	279920 308524	illumigene C. difficile External Control	Positive and Negative Control Reagents for use with the illumigene C. difficile kit	18 tests
Meridian Bioscience	280350 308525	illumigene Group B Streptococcus	DNA Amplification Assay for the Detection of Group B Streptococcus in Vaginal/Rectal Antepartum specimens.	50 tests
Meridian Bioscience	279900 308529	illumigene GBS External Control	Positive and Negative Control Reagents for use with the illumigene GBS kit	18 tests
Meridian Bioscience	280150	illumigene Group A Streptococcus	DNA Amplification Assay for the Detection of Streptococcus pyogenes in throat swab samples. The target DNA sequence (206bp) resides in a region of the pyrogenic exotoxin B (speB).	50 tests
Meridian Bioscience	279910	illumigene GAS External Control	Positive and Negative Control Reagents for use with the illumigene Group A Streptococcus kit	18 tests
Meridian Bioscience	280550	illumigene Mycoplasma	DNA Amplification Assay for the Detection of Mycoplasma pneumoniae in human throat and nasopharyngeal swab specimens.	50 tests
Meridian Bioscience	280250 308540	illumigene Mycoplasma Direct	DNA Amplification Assay for the Detection of Mycoplasma pneumoniae in human throat swabs specimens. Direct detection without DNA extraction.	50 tests
Meridian Bioscience	279940 308541	illumigene Mycoplasma External Control	External Control Materials for use with illumigene Mycoplasma DNA Amplification Assay.	18 tests
Meridian Bioscience	279890 308541	illumigene Mycoplasma Direct External Control	External Control Materials for use with illumigene Mycoplasma Direct DNA Amplification Assay.	18 tests
Meridian Bioscience	280550	illumigene Pertussis	DNA Amplification Assay for the Detection of Bordetella pertussis in nasopharyngeal swab samples	50 tests
Meridian Bioscience	279940	illumigene Pertussis External Control	External Control Materials for use with illumigene Pertussis DNA Amplification Assay.	18 tests
Meridian Bioscience	280450 308545	illumigene Chlamydia (CT)	DNA Amplification Assay for the detection of Chlamydia trachomatis in Endocervical swabs, self-collected vaginal swabs, male urine and female urine	50 tests
Meridian Bioscience	280850 308550	illumigene Gonorrhea (NG)	DNA Amplification Assay for the detection of Neisseria gonorrhoeae in Endocervical swabs, self-collected vaginal swabs, male urine and female urine	50 tests
Meridian Bioscience	281110 308551	M-prep DNA Sample Preparation	DNA sample preparation kit for use with illumigene DNA Amplification Assays	50 units
Meridian Bioscience	6E028N.MER 308552	Mswab	Collection, Preservation and transport system for molecular and culture applications	50 units
Meridian Bioscience	280650 308555	illumigene HSV 1&2	DNA Amplification Assay for the detection and differentiation of herpes simplex virus type 1 (HSV-1) and herpes simplex virus type 2 (HSV-2) in cutaneous and mucocutaneous lesion specimens. 25 HSV-1 + 25 HSV-2 tests	25+25 tests

illumigene System

<i>Meridian Bioscience</i>	279960 308556	illumigene HSV 1&2 External Controls	External Control Materials for use with illumigene HSV 1&2 DNA Amplification Assay.	36+36 tests
<i>Meridian Bioscience</i>	280925 308560	illumigene Malaria	DNA amplification assay for the direct detection of Plasmodium spp. DNA in human whole blood samples	25 tests
<i>Meridian Bioscience</i>	281125 308561	illumigene Malaria PLUS	DNA amplification assay for the direct detection of Plasmodium spp. DNA in human whole blood samples. Including M-Prep extraction columns	25 tests
<i>Meridian Bioscience</i>	279970 308562	illumigene Malaria External Control	Positive and negative control reagents for use with the illumigene Malaria assays	18 tests
<i>Meridian Bioscience</i>	610172 308501	illumipro-10	Automated isothermal amplification and detection system	1 pce

Bacteriology

Bordetella pertussis

ELISA - Antibody Detection

Virotech	EC115.00 341084	Bordetella pertussis IgG/IgA	96 tests
Virotech	EC215G00 341130	Pertussis Toxin IgG	96 tests
Virotech	EC215M00 341132	Pertussis Toxin IgM	96 tests
Virotech	EC215A00 341134	Pertussis Toxin IgA	96 tests

Line Immunoblot

Virotech	WE116G32 343040	Bordetella pertussis + CatACT LINE IgG	32 tests
Virotech	WE116G96 343042	Bordetella pertussis + CatACT LINE IgG	96 tests
Virotech	WE116A32 343046	Bordetella pertussis + CatACT LINE IgA	32 tests
Virotech	WE116A96 343048	Bordetella pertussis + CatACT LINE IgA	96 tests

Borrelia

ELISA - Antibody Detection

Virotech	EC022G00 341002	Borrelia afzelii + VlsE IgG	Antigen mix of B. afzelii strain PKo, B. garinii strain PBr, B. burgdorferi strain ZS7 and VISE	96 tests
Virotech	EC022M00 341006	Borrelia afzelii IgM	Antigen mix of B. afzelii strain PKo, B. garinii strain PBr and B. burgdorferi strain ZS7	96 tests
Virotech	EC023G00 341010	Borrelia VlsE IgG	VlsE of B. burgdorferi s.s. and B. garinii	96 tests
Virotech	EC024G00 341012	Borrelia afzelii + VlsE Europe IgG	Antigen mix of B. afzelii strain PKo, B. garinii strain PBr, B. burgdorferi strain ZS7 and VISE. Special cut-off adjustment	96 tests
Virotech	DC122.00 343100	Borrelia burgdorferi Veterinary ELISA	Semiquantitative and qualitative detection of Borrelia burgdorferi (Species B. burgdorferi sensu stricto) specific IgG-antibodies in dogs- or horses sera	96 tests
Virotech	DC122K80 343101	Borrelia burgdorferi Veterinary IgM Set (anti-dog)	IgM conjugate and IgM cut-off control	1 set

Line Immunoblot

Virotech	WE222G32 343080	Borrelia LINE IgG	VlsE-Mix, p39 (BmpA), p83, iv1 (BBA36), iv2 (BB0323), iv3 (Crasp3), iv4 (pg)	32 tests
Virotech	WE222G96 343081	Borrelia LINE IgG		96 tests
Virotech	WE222M32 343082	Borrelia LINE IgM	OspC, VlsE-Mix, p39 (BmpA), EBV	32 tests
Virotech	WE222M96 343084	Borrelia LINE IgM		96 tests
Virotech	WE223G32 343086	Borrelia in vivo LINE plus IgG	VlsE-Mix, p39 (BmpA), p83, iv1 (BBA36), iv2 (BB0323), iv3 (Crasp3), iv4 (pg), TpN17	32 tests

Bacteriology

Borrelia					
<i>Virotech</i>	WE223G96 343088	Borrelia in vivo LINE plus IgG			96 tests
<i>Virotech</i>	WE224G32 343090	Borrelia Europe LINE IgG	OspC, VlsE-Mix, p39, DbpA-Mix, DbpA-Pko, p58, p83		32 tests
<i>Virotech</i>	WE224G96 343092	Borrelia Europe LINE IgG			96 tests
<i>Virotech</i>	WE224M32 343094	Borrelia Europe LINE IgM	OspC, VlsE-Mix, p39, DbpA-Mix, EBV		32 tests
<i>Virotech</i>	WE224M96 343096	Borrelia Europe LINE IgM			96 tests
<i>Virotech</i>	WE225G32 343098	Borrelia Europe plus TpN17 LINE IgG	OspC, VlsE-Mix, p39, DbpA-Mix, DbpA-Pko, p58, p83, TpN17		32 tests
<i>Virotech</i>	WE225G96 343099	Borrelia Europe plus TpN17 LINE IgG			96 tests
<i>Virotech</i>	DE226G32 343102	Borrelia Veterinary plus OspA LINE, IgG	For dogs and horses. VlsE-Mix DOG (Bb (B31); Bg (IP90)), OspA-Mix (Ba (Pko); Bb (ZS7); Bg (ZQ1)), DbpA-Mix (Ba (Pko); Bg (PBr); Bbav (Pbi)), OspC-Mix (Ba (Pko); Bb (ZS7); Bbav (Pbi)), p39/BmpA (Ba (Pko)), p58 (Bbav (Pbi)), p83 (Ba (Pko)), VlsE-Mix HORSE (Bb (B31); Bg (IP90))		32 tests
<i>Virotech</i>	DE226G96 343104	Borrelia Veterinary plus OspA LINE, IgG			96 tests
<i>Virotech</i>	DE226K62 343106	Borrelia Veterinary plus OspA LINE IgG, Horse Set			1 set
Brucella					
ELISA - Antibody Detection					
<i>Virotech</i>	EC101.00 341014	Brucella IgG/IgM			96 tests
<i>Virotech</i>	EC101.08 341016	Brucella IgA Set			96 tests
Campylobacter					
ELISA - Antigen Detection					
<i>Meridian Bioscience</i>	618096	Premier Campy	Detection of Campylobacter antigens (C. jejuni, C. Coli) in stool specimens		96 tests
Rapid Test - Antigen Detection					
<i>Meridian Bioscience</i>	751530	ImmunoCard STAT! Campy	One step Immuno-Chromatographic Assay for the detection of Campylobacter antigens (C. jejuni, C. coli) in stool specimens		30 tests
Chlamydia					
ELISA - Antibody Detection					
<i>NovaTec</i>	CHLA0070	Chlamydia trachomatis IgA	Highly purified C. trachomatis LGV Type II antigens	strain 434	96 tests
<i>NovaTec</i>	CHLG0070	Chlamydia trachomatis IgG	Highly purified C. trachomatis LGV Type II antigens	strain 434	96 tests
<i>NovaTec</i>	CHLM0070	Chlamydia trachomatis IgM	Highly purified C. trachomatis LGV Type II antigens	strain 434	96 tests

Bacteriology

Chlamydia

NovaTec	CHLA0510	Chlamydia pneumoniae IgA	Highly purified C. pneumoniae CWL-029 antigens	96 tests
NovaTec	CHLG0510	Chlamydia pneumoniae IgG	Highly purified C. pneumoniae CWL-029 antigens	96 tests
NovaTec	CHLM0510	Chlamydia pneumoniae IgM	Highly purified C. pneumoniae CWL-029 antigens	96 tests

Clostridium difficile

ELISA - Antigen Detection

Meridian Bioscience	611096 303100	Premier C. difficile GDH	Detection of C. difficile Common Antigen in stool specimens	96 tests
Meridian Bioscience	616096 303101	Premier Toxins A&B	Microtiter EIA for the detection of C. diff. Toxins A&B in stool specimens	96 tests

Rapid Test - Antigen Detection

Meridian Bioscience	716050 303150	ImmunoCard C. difficile GDH	Rapid enzyme immunoassay for the detection of Clostridium difficile glutamate dehydrogenase (GDH) in fecal specimens	50 tests
Meridian Bioscience	712050 303151	ImmunoCard Toxins A&B	Rapid membrane EIA for the detection of C. diff. Toxins A&B	50 tests

Clostridium tetani

ELISA - Antibody Detection

Progen Biotechnik	7702010 202003	IMMUNOZYM Tetanus	Quantitative determination of anti-tetanus toxoid antibodies in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
Virotech	EC124.00 341112	Tetanus IgG		96 tests

Corynebacterium diphtheriae

ELISA - Antibody Detection

Progen Biotechnik	7702040 202004	IMMUNOZYM Diphtheria	Quantitative determination of anti-diphtheria toxoid antibodies in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
Virotech	EC129.00 341114	Diphtheria IgG		96 tests

Coxiella burnetii

ELISA - Antibody Detection

NovaTec	COX1G0600	Coxiella Burnetii (Q-Fever) Phase 1 IgG	Qualitative determination of IgG antibodies against Coxiella burnetii phase 1 in the early stages of infection in human serum	96 tests
NovaTec	COX2G0600	Coxiella Burnetii (Q-Fever) Phase 2 IgG	Qualitative determination of IgG antibodies against Coxiella burnetii phase 2 in human serum	96 tests
NovaTec	COX2M0600	Coxiella Burnetii (Q-Fever) Phase 2 IgM	Qualitative determination of IgM antibodies against Coxiella burnetii phase 2 in human serum	96 tests

Bacteriology

Entero Hemorrhagic Escherichia coli

ELISA - Antigen Detection

<i>Meridian Bioscience</i>	608096 303201	Premier EHEC	Microtiter EIA for the detection of Shiga-like Toxins I & II in stool and food specimens	96 tests
----------------------------	------------------	---------------------	--	----------

Rapid Test - Antigen Detection

<i>Meridian Bioscience</i>	751630 303203	ImmunoCard STAT! EHEC	Rapid immunoassay for the detection of E. coli Shiga Toxins 1&2 (with differentiation) from enriched stool specimens	30 tests
----------------------------	------------------	------------------------------	--	----------

Reagents (Meridian)

<i>Meridian Bioscience</i>	7516MC2 303204	McConkey Broth	For use with ImmunoCard STAT! EHEC for enrichment of stool specimens	20 x 5ml
----------------------------	-------------------	-----------------------	--	----------

Haemophilus influenza type B (HiB)

ELISA - Antibody Detection

<i>Progen Biotechnik</i>	7704010 202005	IMMUNOZYM HiB IgG	Quantitative determination of IgG antibodies against Haemophilus influenzae type B in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
--------------------------	-------------------	--------------------------	--	----------

Helicobacter pylori

ELISA - Antibody Detection

<i>Virotech</i>	EC143G00 341118	Helicobacter pylori IgG		96 tests
<i>Virotech</i>	EC143A00 341120	Helicobacter pylori IgA		96 tests

ELISA - Antigen Detection

<i>Meridian Bioscience</i>	601396 303002	Premier Platinum HpSA PLUS	Monoclonal based microtiter EIA for the detection of H. pylori Specific Antigens in stool specimens	96 tests
----------------------------	------------------	-----------------------------------	---	----------

Rapid Test - Antigen Detection

<i>Meridian Bioscience</i>	750705 303050	ImmunoCard STAT! HpSA HD	Second generation rapid immunoassay for the detection of H. pylori specific antigens in stool specimens	5 tests
<i>Meridian Bioscience</i>	750720 303051	ImmunoCard STAT! HpSA HD	Second generation rapid immunoassay for the detection of H. pylori specific antigens in stool specimens	20 tests

Rapid Test - Antibody Detection

<i>Meridian Bioscience</i>	710030 303055	ImmunoCard H. pylori	Rapid enzyme immunoassay for IgG detection in serum and whole blood	30 tests
<i>Quidel</i>	0W009 211020	QuickVue H. pylori		10 tests

Line Immunoblot

<i>Virotech</i>	WE243A32 343028	Helicobacter pylori LINE IgA		32 tests
<i>Virotech</i>	WE243A96 343030	Helicobacter pylori LINE IgA		96 tests
<i>Virotech</i>	WE243G32 343034	Helicobacter pylori LINE IgG		32 tests
<i>Virotech</i>	WE243G96 343036	Helicobacter pylori LINE IgG		96 tests

Bacteriology

Legionella

ELISA - Antibody Detection

Meridian Bioscience	633096	Premier L. pneumophila	EIA for the detection of antibodies against L. pneumophila (serogr. 1-6) in human serum	96 tests
IFA - Antibody Detection				
Meridian Bioscience	4873171 303411	L. pneumophila serogroups 1-8	Complete kit with 20 slides (wells 1-5 are coated with inactivated bacteria of serogroup 1, 4, 6, 8; wells 6-10 with inact. bacteria of serogroup 2, 3, 5, 7). The antigens are from european origin.	100 tests
Meridian Bioscience	4873131	L. pneumophila serogroups 1-6	Complete kit with 15 slides with 10 antigen wells each: Antigen wells are coated with inactivated bacteria of L. pneumophila serogroups 1-6. The antigens are from european origin.	150 tests
Meridian Bioscience	4884161	L. pneumophila polyvalent Serogroups 1-6		10 tests
Meridian Bioscience	4873891	L. pneumophila polyvalent Serogroups 1-8	(Pool I (1, 4, 6, 8), Pool II (2, 3, 5, 7))	10 tests
Meridian Bioscience	4873001	L. pneumophila serogroup 1	1 slide with 10 wells containing serogroup 1 antigen	10 tests
Meridian Bioscience	4873002	L. pneumophila serogroup 2	1 slide with 10 wells containing serogroup 2 antigen	10 tests
Meridian Bioscience	4873003	L. pneumophila serogroup 3	1 slide with 10 wells containing serogroup 3 antigen	10 tests
Meridian Bioscience	4873004	L. pneumophila serogroup 4	1 slide with 10 wells containing serogroup 4 antigen	10 tests
Meridian Bioscience	4873005	L. pneumophila serogroup 5	1 slide with 10 wells containing serogroup 5 antigen	10 tests
Meridian Bioscience	4873006	L. pneumophila serogroup 6	1 slide with 10 wells containing serogroup 6 antigen	10 tests
Meridian Bioscience	4873141	L.pneumophila negative control serum		1 ml
Meridian Bioscience	4873011	L.pneumophila positive control serum		1 ml

Rapid Test - Antigen Detection

Meridian Bioscience	751930	TRU Legionella	Rapid immunoassay for the detection of Legionella pneumophila serogroup 1 antigen in human urine specimens	32 tests
---------------------	--------	-----------------------	--	----------

Mycoplasma pneumoniae

ELISA - Antibody Detection

Virotech	EC114.00 341080	Mycoplasma pneumoniae IgG/IgM		96 tests
Virotech	EC114.08 341082	Mycoplasma pneumoniae IgA-Set		96 tests
Virotech	EC214.00 341128	Mycoplasma pn. (rec.) IgG/IgA+M		96 tests

Rapid Test - Antibody Detection

Bacteriology

Mycoplasma pneumoniae

<i>Meridian Bioscience</i>	709030 303551	ImmunoCard Mycoplasma	Rapid enzyme immunoassay for M. pneumoniae IgM detection in human serum	30 tests
Line Immunoblot				
<i>Virotech</i>	WE214G16 343068	Mycoplasma pneumoniae LINE IgG		16 tests
<i>Virotech</i>	WE214M16 343072	Mycoplasma pneumoniae LINE IgM		16 tests
<i>Virotech</i>	WE214A16 343076	Mycoplasma pneumoniae LINE IgA		16 tests

Streptococcus

Rapid Test - Antigen Detection

<i>Meridian Bioscience</i>	752330	TRU Strep Pneumo	A Rapid Immunoassay for the Detection of Streptococcus pneumoniae in Human Urine and Cerebrospinal Fluid Specimens	32 tests
<i>Teomed</i>	241020	Teomed Strep A Test	Rapid detection of Group A Streptococcal antigen from patient throat swab specimens using flocced swabs	20 tests
<i>Quidel</i>	00343IN 211040	QuickVue In-Line Strep A		25 tests
<i>Quidel</i>	20122IN 211045	QuickVue+ Strep A		25 tests

Syphilis

Line Immunoblot

<i>Virotech</i>	WE150G16 343056	Treponema pallidum LINE IgG	Line Immunoblot for the qualitative detection of IgG antibodies to Treponema pallidum antigens TpN47, TmpA (TpN44,5), TpN17, TpN15 in human serum. These markers are highly specific for primary, secondary and latent syphilis.	16 tests
<i>Virotech</i>	WE150G32 343058	Treponema pallidum LINE IgG		32 tests
<i>Virotech</i>	WE150M16 343062	Treponema pallidum LINE IgM	Line Immunoblot for the qualitative detection of IgM antibodies to Treponema pallidum antigens TpN47, TmpA (TpN44,5), TpN17, TpN15 in human serum. These markers are highly specific for primary, secondary and latent syphilis.	16 tests
<i>Virotech</i>	WE150M32 343064	Treponema pallidum LINE IgM		32 tests

Yersinia enterocolitica

ELISA - Antibody Detection

<i>Virotech</i>	EC142.00 341116	Yersinia enterocolitica IgG/IgA		96 tests
-----------------	--------------------	--	--	----------

Agglutination Tests

<i>Progen Biotechnik</i>	PR87025 202508	Yersinia enterocolitica 0:3	Monoclonal antibody for identification of bacteria colonies of Yersinia enterocolitica serogroup 0:3	40 tests
<i>Progen Biotechnik</i>	PR87035 202509	Yersinia enterocolitica 0:9	Monoclonal antibody for identification of bacteria colonies of Yersinia enterocolitica serogroup 0:9	40 tests

Bacteriology

Yersinia enterocolitica				
Progen Biotechnik	PR87015	Yersinia enterocolitica O:3/O:9	Monoclonal antibody agglutination test for identification of bacterial colonies of Yersinia enterocolitica serotypes O:3 and O:9	40 tests
Line Immunoblot				
Virotech	WE242A32 343012	Yersinia enterocolitica LINE IgA	Qualitative detection of specific antibodies directed against antigens of the 70kb virulence plasmid of pathogenic Yersinia species. Includes recombinant antigens of Yersinia outer membrane proteins Yop M, Yop H, Yop B, V-Ag, Yop D, Yop E and Yop DD	32 tests
Virotech	WE242A96 343014	Yersinia enterocolitica LINE IgA		96 tests
Virotech	WE242G32 343018	Yersinia enterocolitica LINE IgG		32 tests
Virotech	WE242G96 343020	Yersinia enterocolitica LINE IgG		96 tests
Virotech	WE242M32 343024	Yersinia enterocolitica LINE IgM		32 tests

Virology

Adenovirus					
ELISA - Antibody Detection					
<i>Virotech</i>	EC121.00 341108	Adenovirus IgG/IgM			96 tests
<i>Virotech</i>	EC121.08 341110	Adenovirus IgA-Set			96 tests
ELISA - Antigen Detection					
<i>Meridian Bioscience</i>	696006	Premier Adenoclone 40/41	Microtiter EIA for direct detection of Adenoviruses serotypes 40/41 in stool samples		48 tests
<i>Meridian Bioscience</i>	696007	Premier Adenoclone	Microtiter EIA for detection of adenovirus antigen in stool samples and for confirmation in cell culture		48 tests
Rapid Test - Antigen Detection					
<i>Meridian Bioscience</i>	751030 300651	ImmunoCard STAT! Adenovirus	Rapid Immunoassay for the detection of Adenovirus antigens		30 tests
Chikungunya Virus					
ELISA - Antibody Detection					
<i>NovaTec</i>	CHIG0590	Chikungunya IgG	Chikungunya antigens		96 tests
<i>NovaTec</i>	CHIM0590	Chikungunya IgM	μ-capture, Chikungunya antigens		96 tests
CMV					
ELISA - Antibody Detection					
<i>Virotech</i>	EC113.00 341076	Cytomegalo Virus IgG/IgM			96 tests
Dengue Virus					
ELISA - Antibody Detection					
<i>NovaTec</i>	DENG0120	Dengue IgG	Purified Dengue Virus 2 antigens from strain 16681		96 tests
<i>NovaTec</i>	DENM0120	Dengue IgM	Purified Dengue Virus 2 antigens from strain 16681		96 tests
EBV					
ELISA - Antibody Detection					
<i>Virotech</i>	EC102.00 341018	Epstein Barr Virus IgG/IgM	Semiquantitative and qualitative detection of IgM- and IgG-antibodies against Epstein Barr Virus		96 tests
<i>Virotech</i>	EC203M00 341022	EBV VCA IgM	Semiquantitative		96 tests
<i>Virotech</i>	EC204.00 341024	EBV EBNA IgG	Semiquantitative		96 tests
<i>Virotech</i>	EC202.00 341026	EBV Early IgG	Semiquantitative		96 tests
<i>Virotech</i>	EC205.00 341028	EBV VCA IgG	Semiquantitative		96 tests
IFA - Antibody Detection					

Virology

EBV				
Meridian Bioscience	EB-100 300011	Merifluor EBV IgG	Complete kit. Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EB-150 300012	Merifluor EBV IgM	Complete kit. Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EN-100 300013	Merifluor EBNA ACIF	Complete kit. Slides coated with 20-40% EBNA pos. Burkitt lymphoma Raji cell line and with EBNA neg. cells	10x10 tests
Meridian Bioscience	EA-101 300014	Merifluor EBV-EA	Complete kit. Slides coated with Raji cells from Burkitt lymphoma	10x10 tests
Meridian Bioscience	EB-211 300021	EBV IgM/IgG slides pack	Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EA-211 300024	EBV-EA slides pack	Slides coated with Raji cells from Burkitt lymphoma	10x10 tests
Agglutination Tests				
Meridian Bioscience	776060 300051	Monospot	Slide test using fresh horse erythrocytes	20 tests
Meridian Bioscience	776150 300052	Monospot Latex	Latex particle agglutination assay	50 tests
Rapid Test - Antigen Detection				
Quidel	20121 211050	QuickVue+ Infect. Mononucleosis		20 tests
Rapid Test - Antibody Detection				
Meridian Bioscience	751830 300041	TRU EBV-G	One step Immuno-Chromatographic assay for determination of anti-VCA IgG and anti-EBNA-1 IgG	30 tests
Meridian Bioscience	776500	Monolert	Rapid EIA test for anti-EBNA-1 IgG and IgM status	100 tests
Meridian Bioscience	776520	Monolert Control Sera	Positive and negative control sera	25 tests
Line Immunoblot				
Virotech	WE102G32 343002	EBV LINE IgG	Qualitative detection of IgG antibodies to EBNA1, gp125, p18 and EA-D in human serum. High sensitivity and specificity by using synthetical peptides. Clear differentiation between primary and past infection due to the double-cut-off-strategy.	32 tests
Virotech	WE102G96 343004	EBV LINE IgG		96 tests
Virotech	WE102M96 343010	EBV LINE IgM		96 tests
Virotech	WE102M32 343008	EBV LINE IgM	Qualitative detection of IgM antibodies to gp125, p18 and EA-D in human serum. High sensitivity and specificity by using synthetical peptides.	32 tests
Enterovirus				
ELISA - Antibody Detection				
Virotech	EC116A00 341086	Enterovirus IgA		96 tests
Virotech	EC116G00 341088	Enterovirus IgG		96 tests
Virotech	EC116M00 341090	Enterovirus IgM		96 tests

Virology

FSME (TBE; Tick-Borne Encephalitis)

ELISA - Antibody Detection

<i>Virotech</i>	EC117.00 341092	TBE IgG/IgM		96 tests
-----------------	--------------------	--------------------	--	----------

Hanta Virus

ELISA - Antibody Detection

<i>Progen Biotechnik</i>	PR59165	Hantavirus (Dobrava/Hantaan) IgG/IgM	Semiquantitative Elisa for the detection of IgG and IgM antibodies against hantavirus (Dobrava and Hantaan serotypes) in human serum	96 tests
<i>Progen Biotechnik</i>	PR59156	Hantavirus (Puumala) IgG/IgM	Semiquantitative Elisa for the detection of IgG and IgM antibodies against hantavirus (Puumala group) in human serum	96 tests

IFA - Antibody Detection

<i>Progen Biotechnik</i>	PR77065 202504	Hantavirus (Dobrava/Hantaan) Antibody IFA	For the detection of antibodies against Hantavirus Hantaan and Dobrava serotypes in human serum.	60 tests
<i>Progen Biotechnik</i>	PR77056 202505	Hantavirus (Puumala) Antibody IFA	For the detection of antibodies against hantavirus (Puumala group) in human serum	60 tests
<i>Progen Biotechnik</i>	PR77060 202506	Hantavirus (Seoul) Antibody IFA	For the detection of antibodies against hantavirus (Seoul group) in human serum	60 tests
<i>Progen Biotechnik</i>	FR77090 203110	Hantavirus IF Reagent Kit Rat/Mouse	For use in an indirect immunofluorescence test for the the determination of rat or mouse antibodies against hantavirus (hantaan/puumala/seoul) in rat and mouse serum.	60 tests
<i>Progen Biotechnik</i>	PR97065	Slide Pack (Dobrava/Hantaan)		60 tests
<i>Progen Biotechnik</i>	PR97056	Slide pack (Puumala)		60 tests
<i>Progen Biotechnik</i>	PR97060 203111	Slide pack (Seoul)		60 tests

HSV

ELISA - Antibody Detection

<i>Virotech</i>	EC108.00 341030	Herpes simplex Virus Screen IgG/IgM		96 tests
<i>Virotech</i>	EC130.00 341036	Herpes simplex Virus I (gG1) spec. IgG/IgM		96 tests
<i>Virotech</i>	EC131.00 341040	Herpes simplex Virus II (gG2) spec. IgG/IgM		96 tests

Line Immunoblot

<i>Virotech</i>	WE130G16 343052	HSV LINE IgG		16 tests
-----------------	--------------------	---------------------	--	----------

Influenza Virus

ELISA - Antibody Detection

<i>Virotech</i>	EC118.00 341098	Influenza A IgG/IgM		96 tests
-----------------	--------------------	----------------------------	--	----------

Virology

Influenza Virus				
<i>Virotech</i>	EC118A00 341100	Influenza A IgA		96 tests
<i>Virotech</i>	EC119.00 341102	Influenza B IgG/IgM		96 tests
<i>Virotech</i>	EC119A00 341104	Influenza B IgA		96 tests
Rapid Test - Antigen Detection				
<i>Meridian Bioscience</i>	751230	TRU FLU	Rapid Immunoassay with separate A / B results	32 tests
<i>Meridian Bioscience</i>	751110	FLU/RSV Positive Control	3 vials of 8 tests each	24 tests
<i>Quidel</i>	20183IN 211036	QuickVue Influenza A+B		25 tests
Measles (Rubeola) Virus				
ELISA - Antibody Detection				
<i>Virotech</i>	EC105G00 341044	Measles IgG		96 tests
<i>Virotech</i>	EC105M00 341046	Measles IgM		96 tests
Mumps Virus				
ELISA - Antibody Detection				
<i>Virotech</i>	EC106.00 341052	Mumps IgG/IgM		96 tests
Norovirus				
Rapid Test - Antigen Detection				
<i>Meridian Bioscience</i>	757020 300752	ImmunoCard STAT! Norovirus	Rapid Immunoassay for the detection of Norovirus antigens in stool specimens	20 tests
<i>Meridian Bioscience</i>	757001 300759	ImmunoCard STAT! Norovirus positive control		1 test
Parainfluenza Virus				
ELISA - Antibody Detection				
<i>Virotech</i>	EC148.00 341122	ParaScreen IgG/IgM		96 tests
<i>Virotech</i>	EC148.08 341124	ParaScreen IgA-Set		96 tests
Parvovirus				
ELISA - Antibody Detection				
<i>NovaTec</i>	PARG0370 321511	Parvovirus B19 IgG	recombinant Parvovirus B19 antigens	96 tests
<i>NovaTec</i>	PARM0370 321512	Parvovirus B19 IgM	recombinant Parvovirus B19 antigens	96 tests

Virology

Rotavirus

ELISA - Antigen Detection

<i>Meridian Bioscience</i>	696004	Premier Rotaclone	Using Rotavirus monoclonal antibodies	48 tests
----------------------------	--------	--------------------------	---------------------------------------	----------

Rapid Test - Antigen Detection

<i>Meridian Bioscience</i>	750030 300551	ImmunoCard STAT! Rotavirus	Rapid Immunoassay for the detection of Rotavirus antigen	30 tests
<i>Meridian Bioscience</i>	751120 300552	Rapid Strip Rota / Adeno	Rapid Immunoassay with separate bands for Rotavirus and Adenovirus	20 tests

RSV

ELISA - Antibody Detection

<i>Virotech</i>	EC107.00 341056	Respiratory syncytial Virus IgG/IgM		96 tests
<i>Virotech</i>	EC107.08 341058	Respiratory syncytial Virus IgA-Set		96 tests

Rapid Test - Antigen Detection

<i>Meridian Bioscience</i>	751110	FLU/RSV Positive Control	3 vials of 8 tests each	24 tests
<i>Meridian Bioscience</i>	751330 300453	TRU RSV	Rapid Immunoassay for the detection of Respiratory Syncytial Virus (RSV) antigens	32 tests
<i>Quidel</i>	20193 211056	QuickVue RSV		20 tests

Rubella Virus

ELISA - Antibody Detection

<i>Virotech</i>	EC109L00 341060	Rubella virus IgG CSF ELISA		96 tests
-----------------	--------------------	------------------------------------	--	----------

VZV

ELISA - Antibody Detection

<i>Virotech</i>	EC110.00 341062	Varicella zoster Virus IgG/IgM		96 tests
<i>Virotech</i>	EC110.08 341064	Varicella zoster Virus IgA-Set		96 tests

IFA - Antigen Detection

<i>Meridian Bioscience</i>	500113 300210	Merifluor VZV	For detection of VZV in specimens and culture	60 tests
----------------------------	------------------	----------------------	---	----------

Mycology

Candida albicans

ELISA - Antibody Detection

<i>Virotech</i>	EC111.00 341072	Candida albicans IgG/IgM		96 tests
<i>Virotech</i>	EC111.08 341074	Candida albicans IgA-Set		96 tests

Cryptococcus neoformans

ELISA - Antigen Detection

<i>Meridian Bioscience</i>	602096	Premier Cryptococcal Antigen	EIA for the detection of Cr. neoformans antigen in serum and CSF	192 tests
----------------------------	--------	-------------------------------------	--	-----------

Agglutination Tests

<i>Meridian Bioscience</i>	140100 308151	Cryptococcal Latex Aggl. System (CALAS)	For the detection of Cr. neoformans Ag in serum and CSF	100 tests
<i>Meridian Bioscience</i>	140050	CALAS Pronase Reagent	Enzyme for the treatment of serum samples to remove aspecific interferences	5 vials

Fungal Serology

ELISA - Antibody Detection

<i>Meridian Bioscience</i>	603096	Premier Coccidioides	Detection of Anti-C. immitis IgG and IgM in serum or CSF	96 tests
----------------------------	--------	-----------------------------	--	----------

Immunodiffusion

<i>Meridian Bioscience</i>	100201 308031	Histoplasma Antigen		1 ml
<i>Meridian Bioscience</i>	100601 308032	Anti-Histoplasma		1 ml
<i>Meridian Bioscience</i>	100401 308034	Coccidioides "F" Antigen		1 ml
<i>Meridian Bioscience</i>	100801 308035	Anti-Coccidioides "F"		1 ml
<i>Meridian Bioscience</i>	103001 308038	Coccidioides "TP" Antigen		1 ml
<i>Meridian Bioscience</i>	103101 308039	Anti-Coccidioides "TP"		1 ml
<i>Meridian Bioscience</i>	101009 308040	Immunodiffusion Plates	4 series/plate	9 plates
<i>Meridian Bioscience</i>	101012	Immunodiffusion Plates	1 series/plate	12 plates

Pneumocystis carinii

IFA - Antigen Detection

<i>Meridian Bioscience</i>	222030 308211	Merifluor Pneumocystis	DFA for detection of Pneumocystis carinii cysts and trophozoites in respiratory specimens	30 tests
----------------------------	------------------	-------------------------------	---	----------

Parasitology

Ascaris

ELISA - Antibody Detection

NovaTec	ASCG0020	Ascaris lumbricoides IgG	PBS extract of Ascaris	96 tests
---------	----------	---------------------------------	------------------------	----------

Cryptosporidium / Giardia / Entamoeba

ELISA - Antigen Detection

NovaTec	GIA0160S	Giardia lamblia Antigen in stool		96 tests
---------	----------	---	--	----------

IFA - Antigen Detection

Meridian Bioscience	250050 306110	Merifluor Crypto & Giardia	Direct FA for Cryptosporidium and Giardia in stool specimens	50 tests
---------------------	------------------	---------------------------------------	--	----------

Rapid Test - Antigen Detection

Meridian Bioscience	750830 306151	ImmunoCard STAT! Crypto / Giardia	Rapid immunoassay for the qualitative detection of Cryptosporidium parvum and Giardia lamblia specific antigens from fresh or preserved (SAF, Cary-Blair etc.) stool specimens	30 tests
---------------------	------------------	--	--	----------

Meridian Bioscience	750820	ImmunoCard STAT! Crypto / Giardia Positive Controls	1 ml Giardia and 1 ml Cryptosporidium Positive Control vials	16 tests
---------------------	--------	--	--	----------

Meridian Bioscience	751420 306153	ImmunoCard STAT! CGE	One-step immunochromatographic test for the qualitative detection of Cryptosporidium parvum, Giardia lamblia and Entamoeba histolytica antigens in human stool samples	20 tests
---------------------	------------------	-----------------------------	--	----------

Meridian Bioscience	751403 306154	ImmunoCard STAT! CGE Positive Control	External Positive Control Set for ImmunoCard STAT! CGE	3 tests
---------------------	------------------	--	--	---------

Meridian Bioscience	757220	Rapid Strip Crypto & Giardia	Rapid Immunoassay for the detection of Cryptosporidium parvum and Giardia lamblia antigens in native stool specimens	20 tests
---------------------	--------	---	--	----------

Leishmania

ELISA - Antibody Detection

NovaTec	LEIG0310	Leishmania IgG	Purified Leishmania antigens	96 tests
---------	----------	-----------------------	------------------------------	----------

NovaTec	LEIG0310V	Leishmania IgG canine	Purified Leishmania antigens	96 tests
---------	-----------	------------------------------	------------------------------	----------

Toxoplasma gondii

ELISA - Antibody Detection

NovaTec	TOXM0460	Toxoplasma gondii IgM	μ-Capture, P30 antigen	96 tests
---------	----------	------------------------------	------------------------	----------

NovaTec	TOXG0460	Toxoplasma gondii IgG	quantitative, Toxoplasma gondii (Lysate) antigens	96 tests
---------	----------	------------------------------	---	----------

NovaTec	TOXGA460	Toxoplasma gondii IgG Avidity Reagents	IgG avidity reagents intended to differentiate between past and acute infection. To be used with TOXG0460 Toxoplasma gondii IgG	96 tests
---------	----------	---	---	----------

Parasitology

Trichinella

ELISA - Antibody Detection

NovaTec	TRIG0480	Trichinella spiralis IgG	Trichinella antigens	96 tests
---------	----------	---------------------------------	----------------------	----------

Trypanosoma Cruzi (Chagas)

ELISA - Antibody Detection

NovaTec	CHAG0560	Chagas (Trypanosoma cruzi) IgG	Recombinant Trypanosoma cruzi antigens (TcF)	96 tests
---------	----------	---------------------------------------	--	----------

Closed Collection & Concentration Systems

Parapak

Collection and Processing Systems for Stool Specimens

Meridian Bioscience	970120	Para-Pak Macro-CON	stool concentration system	120 units
Meridian Bioscience	900312 307001	Para-Pak Clean vial	empty vial	120 units
Meridian Bioscience	900412	Para-Pak 10% Formalin	10% neutral buffered formalin	120 units
Meridian Bioscience	900212 307002	Para-Pak SAF	sodium acetate formalin	120 units
Meridian Bioscience	901312	Para-Pak Ecofix	new safer fixative without formalin	120 units
Meridian Bioscience	900112	Para-Pak LV-PVA	polyvinyl alcohol with mercuric chloride	120 units
Meridian Bioscience	901012	Para-Pak Zn-PVA	mercury free PVA (with zinc sulphate)	120 units
Meridian Bioscience	900612 307003	Para-Pak C&S	transport medium for enteric pathogens	120 units

Collection Systems with Integrated Filters

Parapak Plus+Ultra

Collection and Processing Systems for Stool Specimens

Meridian Bioscience	980212 307021	Para-Pak ULTRA SAF	closed system for 15/50 ml tubes	120 units
Meridian Bioscience	980412	Para-Pak ULTRA 10% Formalin	closed system for 15/50 ml tubes	120 units
Meridian Bioscience	981312	Para-Pak ULTRA Ecofix	closed system for 15/50 ml tubes, safe fixative	120 units
Meridian Bioscience	800400	Conical tubes and caps	recommended for use with Para-Pak ULTRA, 15 ml	500 units

Filtering Centrifuge Tube

SpinCon

Collection and Processing Systems for Stool Specimens

<i>Meridian Bioscience</i>	972200 307031	SpinCon	Filtering tube for easy concentration, no chemical required	200 units
<i>Meridian Bioscience</i>	972500 307032	SpinCon	Filtering tube for easy concentration, no chemical required	500 units

Reagents for Parasitology

Reagents

Collection and Processing Systems for Stool Specimens

<i>Meridian Bioscience</i>	401301	Ecofix bulk	new safer fixative without formalin	500 ml
<i>Meridian Bioscience</i>	FK4318 307044	Ethyl Acetate	recommended for the Ritchie concentration	225 ml
<i>Meridian Bioscience</i>	801300 307043	Hemo-De	safe substitute for Ethyl Acetate and Xylene	400 ml
<i>Meridian Bioscience</i>	400101 307041	Wheatley's Trichrome Stain	ready-to-use reagent for permanent staining	500 ml
<i>Meridian Bioscience</i>	400501 307042	Para-Pak Trichrome Blue Stain	specific, ready-to-use stain for Microsporidia	500 ml
<i>Meridian Bioscience</i>	801400	Ecostain Set	reagents set for staining Ecofix preparations	1 set

Allergy Diagnostics

Immunoblot

<i>Mediwiss Analytic</i>	120210	AllergyScreen INH 16+T	Immunoblot for the determination of 16 different specific IgE to inhalative allergens (Swiss Panel)	12 tests
<i>Mediwiss Analytic</i>	120220	AllergyScreen NM 16+T	Immunoblot for the determination of 16 different specific IgE to alimentary allergens (Swiss Panel)	12 tests
<i>Mediwiss Analytic</i>	120230	AllergyScreen Panel INH 10	Immunoblot for the determination of 10 different specific IgE to inhalative allergens	12 tests
<i>Mediwiss Analytic</i>	120400	AllergyScreen Panel NM/INH 10	Immunoblot for the determination of specific IgE to 7 alimentary and 3 inhalative allergens	12 tests
<i>Mediwiss Analytic</i>	122005	AllergyScreen Panel 1	Immunoblot for the determination of 20 different specific IgE to inhalative and alimentary allergens	12 tests
<i>Mediwiss Analytic</i>	122010	AllergyScreen Panel 2	Immunoblot for the determination of 20 different specific IgE to inhalative allergens	12 tests
<i>Mediwiss Analytic</i>	122020	AllergyScreen Panel 3	Immunoblot for the determination of 20 different specific IgE to alimentary allergens	12 tests
<i>Mediwiss Analytic</i>	122030	AllergyScreen Panel 4	Immunoblot for the determination of 20 different specific IgE to inhalative and alimentary allergens. Designed for pediatric patients.	12 tests
<i>Mediwiss Analytic</i>	122040	AllergyScreen Panel 5	Immunoblot for the determination of 20 different specific IgE to inhalative allergens and bee/wasp venom.	12 tests
<i>Mediwiss Analytic</i>	A992	CCD-Blocking-Solution	Blocking of cross-reactive carbohydrate determinants (CCDs) for use with AllergyScreen test kits. To be mixed with serum prior to incubation (10% CCD solution + 90% serum)	500 µl
<i>Mediwiss Analytic</i>	120401	Improvio C Reader	Scanner and software for the analysis of the AllergyScreen tests	1 pce
<i>Mediwiss Analytic</i>	120409	Calibration Card	Used to calibrate Improvio scanner	1 pce

Clinical Chemistry

Gastrointestinal Diseases

ELISA - Antigen Detection

<i>Epitope Diagnostics</i>	KT-850	Fecal Human Hemoglobin		96 tests
<i>Epitope Diagnostics</i>	KT849	Fecal Calprotectin qn Elisa	Elisa for the quantitative measurement of human Calprotectin level in stool. Including standards up to 2000 µg/g.	96 tests
<i>Epitope Diagnostics</i>	KT843	Calprotectin Sample Collection Kit	Specially designed for easy collection of a defined amount of fecal sample. Pre-filled with sample extraction buffer.	50 units

Rapid Test - Antigen Detection

<i>Epitope Diagnostics</i>	KT313 230010	EpiTub Fecal Occult Blood Test	Monoclonal based immunochromatographic assay for the detection of fecal occult blood	30 tests
--------------------------------	-----------------	---------------------------------------	--	----------

hCG

Rapid Test - Antigen Detection

<i>Quidel</i>	00178 211015	QuickVue+ One-Step hCG Combo		30 tests
---------------	-----------------	-------------------------------------	--	----------

Lipid Metabolism

ELISA - Antigen Detection

<i>Progen Biotechnik</i>	7700040	IMMUNOZYM Lp(a)	Quantitative Elisa to determine the concentration of Lp(a)	96 tests
------------------------------	---------	------------------------	--	----------

Reagents

Reagents				
Reagents (Meridian)				
Meridian Bioscience	XX-715 309001	GullSORB	IgG inactivation reagent, ready to use	5 ml
Meridian Bioscience	4870141 309002	Mounting medium		3 ml
Meridian Bioscience	4874001 309003	Specimen diluent	Lyophilized (20 ml diluted reagent), recommended for EBV-EA and L. pneumophila IFA	1 vial
Meridian Bioscience	4873691 309004	PBS buffer	Concentrate (30x for 1L working solution)	33 ml
Meridian Bioscience	4873151	FITC-labelled Anti-Human Immunoglobulin	Lyophilized (for 2,5 ml)	1 vial
Reagents (Virotech)				
Virotech	EC250.00 342002	Validation kit	To check ELISA processors	96 tests
Virotech	EN250K60 342004	Pipetting Control-Set	To check the pipetting process in combination with the validation kit	2.5 ml
Virotech	161101 342010	RF-Sorbotech Pre-adsorption	1 vial (2 ml) for 40 samples in human serum or CSF	2 ml
Virotech	161102 342012	RF-Sorbotech Pre-adsorption	1 vial (10 ml) for 200 samples in human serum or CSF	10 ml
Virotech	B/300.11 342014	RF-Sorbotech Pre-adsorption	1 vial RF (2ml), 40 tubes with 450 µl dilution dilution buffer each. For 40 samples in human serum	1 set
Virotech	B/300.00 342016	RF-Sorbotech Pre-adsorption combi pack	2 vials (2ml each) and dilution buffer for 80 samples in human serum or CSF	1 set
Virotech	131101 349001	IgG Conjugate for ELISA (anti-human)		11 ml
Virotech	131201 349002	IgM Conjugate for ELISA (anti-human)		11 ml
Virotech	131301 349003	IgA Conjugate for ELISA (anti-human)		11 ml
Virotech	141101 349004	Citrate-Stopping Solution for ELISA		6 ml
Virotech	141201 349006	PBS-wash buffer (conc 20x) for ELISA		50 ml
Virotech	141301 349007	PBS dilution buffer for ELISA (blue, ready to use)		50 ml
Virotech	151101 349008	TMB substrate solution for ELISA		11 ml
Virotech	141401 349009	VZV IgM dilution buffer (green, ready to use) for ELISA		50 ml
Virotech	ENXXXK40	ELISA IgA controls	IgA negative, cut-off, and positive controls. Only for B. pertussis, Enterovirus, H. pylori, Influenza A/B, Pertussis Toxin and Y. enterocolitica	1 set
Virotech	ENXXXK60	ELISA IgG controls	IgG negative, cut-off, and positive controls	1 set

Reagents

Reagents				
Virotech	ENXXXK80	ELISA IgM controls	IgM negative, cut-off, and positive controls	1 set
Virotech	ECXXXL40	CSF IgA Standards	Available for HSV Screen, VZV. 4 ready to use CSF-Standards for a standard curve.	4 x 1 ml
Virotech	ECXXXL60	CSF IgG Standards	Available for Borrelia, CMV, EBV, HSV Screen, HSV-1, HSV-2, Influenza B, Measles, Mumps, Rubella, TBE, VZV. 4 ready to use CSF-Standards for a standard curve.	4 x 1 ml
Virotech	ECXXXL80	CSF IgM Standards	Available for Borrelia, Measles, TBE, VZV. 4 ready to use CSF-Standards for a standard curve.	4 x 1 ml
Virotech	ENXXXL65	Antibody-Index (AI) Control	Available for Borrelia afzelii+VISE IgG, EBV IgG, VZV IgG, Measles IgG, TBE IgG, Rubella IgG, HSV Screen IgG/IgA	4 x 1 ml
Virotech	WEXXXPXX	LINE positive control		0.5 ml
Virotech	WEXXXNXX	LINE negative control		0.5 ml
Virotech	WEXXXC40	LINE IgA cut-off control (conc.)		0.5 ml
Virotech	WEXXXC60	LINE IgG cut-off control (conc.)		0.5 ml
Virotech	WEXXXC80	LINE IgM cut-off control (conc.)		0.5 ml
Virotech	WE200.08 349102	Serum dilution/wash buffer for LINE (conc. 10x)		50 ml
Virotech	WE200.02 349104	Chromogene/Substrate solution, BCIP/NBT for LINE		57 ml
Virotech	WE400.42 349106	IgA Conjugate for LINE (anti-human, 100x conc.)		0.7 ml
Virotech	WE400.62 349108	IgG Conjugate for LINE (anti-human, 100x conc.)		0.7 ml
Virotech	WE400.82 349110	IgM Conjugate for LINE (anti-human, 100x conc.)		0.7 ml
Virotech	WE300.08 349112	Incubation tray with 8 cavities for LINE		1 pce

Index

Adenovirus	16	Leukemia	6
Allergy Diagnostics	25	Lipid Metabolism	26
Ascaris	22	Luminex Accessories	5
Autoimmune Diagnostics	2	Luminex System	1
Bordetella pertussis	9	Measles (Rubeola) Virus	19
Borrelia	1	Measles, Mumps, Rubella, Varizella	3
	9	Mumps Virus	19
	10	Mycoplasma pneumoniae	13
Brucella	10		14
C9orf72 Gene	6	NGS (Next Generation Sequencing)	6
Campylobacter	10	Norovirus	19
Candida albicans	21	Parainfluenza Virus	19
Chikungunya Virus	16	Parapak	23
Chlamydia	10	Parapak Plus+Ultra	23
	11	Parvovirus	19
Clostridium difficile	11	Pneumocystis carinii	21
Clostridium tetani	11	Quality Control Sera	4
CMV	16	Reagents	24
Corynebacterium diphtheriae	11		27
Coxiella burnetii	11		28
Cryptococcus neoformans	21	Rotavirus	20
Cryptosporidium / Giardia / Entamoeba	22	RSV	20
Dengue Virus	16	Rubella Virus	20
EBV	3	SpinCon	24
	16	Streptococcus	14
	17	Syphilis	14
Entero Hemorrhagic	12	Toxoplasma gondii	22
Escherichia coli		Toxoplasma, Rubella, CMV, HSV	3
Enterovirus	17	Trichinella	23
Fragile X	6	Trypanosoma Cruzi (Chagas)	23
FSME (TBE; Tick-Borne Encephalitis)	18	VZV	20
Fungal Serology	21	Yersinia enterocolitica	14
Gastrointestinal Diseases	26		15
Haemophilus influenza type B (HiB)	12		
Hanta Virus	18		
hCG	26		
Helicobacter pylori	12		
HPV	1		
HSV	3		
	18		
illumigene	7		
	8		
Influenza Virus	18		
	19		
Legionella	13		
Leishmania	22		